Laser Safety: Safe clinical application of laser radiation

Frits F.M. de Mul MEDPHYS Software & Services

2012


www.medphys.nl

email: info(at)medphys.nl

- NOHD: safe distance to avoid eye damage
- Standard distance for safety: 0.25 m
 for 10 s exposure time.
- If laser intensity [W/m²] at the standard distance is too high: protection glasses compulsory.

Damage of excessive laser exposure to eye and skin

	eye			skin	
Wave- length	Photo- chemical	Retina	Thermal	Erythema formation	Thermal
< 400 nm	+	-	+	+	+
400-600	+	+	+	-	+
600-700	-	+	+	-	+
> 700	-	-	+	-	+



- Safety condition at target:
- Intensity $I \leq Maximum Permissible Exposure <math>MPE$
- **I** and **MPE** in W/m²
- **NOHD** is minimum distance where $I \leq MPE$

MPE: maximum permissible exposure


MPE is a function of:


- Wavelength
- Pulse duration (or continuous)
- Pulse frequency
- Angular subtense (viewing angle of light source from target)

MPE-formulae are tabulated in

- Directive 2006/25/EU of European Parliament and Commission
- IEC 60825-1
- NEN 60825-1


MPE: maximum permissible exposure : eye





 $\alpha =$ angular
subtense₆


MPE: maximum permissible exposure : eye




- Safety condition: $I \leq MPE$ in W/m²
- With **I** = **P.PC** / **S**
- $I = intensity [W/m^2]$
- **P** = laser power [W]
- **PC** = laser Profile Correlation Factor [-]
- S = target area [m²]

PC: Profile Correlation Factor

PC depends on laser mode structure




Cylindrical transverse modes



Rectangular transverse modes

Local intensity in modes is larger than averaged over the mode envelope.

→ Profile Correlation Factor **PC**.



Effective target area:

$$S = \frac{\pi}{4} D'^2, \ D' = D + w$$

D' = effective diameter


$$\frac{4S}{\pi} = D'^2 = (D+w)^2$$

Safety condition:

$$I = \frac{P.PC}{S} \le MPE \qquad \Rightarrow \quad S \ge \frac{P.PC}{MPE}$$

$$\Rightarrow S \ge \frac{P.PC}{MPE}$$

$$D = \sqrt{\frac{4S}{\pi}} - w \ge \sqrt{\frac{4P.PC}{\pi.MPE}} - w$$




Effective target area:

$$S = \frac{\pi}{4} D'^2, \quad D' = D + w$$

Safety condition:
$$D \ge \sqrt{\frac{4S}{\pi} - w} = \sqrt{\frac{4P.PC}{\pi.MPE}} - w$$

$$NOHD = R + a = \frac{D}{2 \cdot \tan(\theta/2)} + a =$$

$$= \frac{1}{2 \cdot \tan(\theta/2)} \left[\sqrt{\frac{4P \cdot PC}{\pi \cdot MPE}} - w \right] + a$$




NOHD depends on:

- C
- ≈ 1/θ
- √P
- √PC
- √1/MPE
- -W
- a

$$NOHD = R + a = \frac{D}{2 \cdot \tan(\theta/2)} + a =$$

$$= \frac{1}{2 \cdot \tan(\theta/2)} \left[\sqrt{\frac{4P \cdot PC}{\pi \cdot MPE}} - w \right] + a$$

Focus as seen from the target



NOHD is also dependent on α : "angular subtense"

"angular subtense" = angle to view the source (= focus) from the target position


$$\tan\frac{1}{2}\alpha = \frac{\frac{1}{2}w}{R} = \frac{w}{2(NOHD - a)}$$

A smaller waist will lead to a sharper (intenser) spot on the target

For certain situations (e.g. visible light):


MPE depends on α (tabulated in EU-Directive), and thus on *NOHD*

Two relations between NOHD and MPE:



$$\tan \frac{1}{2}\alpha = \frac{w}{2R} = \frac{w}{2(NOHD - a)}$$

So, NOHD =
$$f(\alpha)$$




But also:

 $NOHD = f(MPE) = f(\alpha)!$


From Safety Condition: *I* ≤ *MPE*

$$NOHD = R + a = \frac{1}{2 \cdot \tan(\theta/2)} \left[\sqrt{\frac{4P \cdot PC}{\pi \cdot MPE}} - w \right] + a$$

Thus: iterative approach necessary !!



Transmission of safety goggles



In case NOHD > safety distance b (from focus): safety goggles necessary!


Transmission T of safety goggles and Optical density OD:

$$T = \left[\frac{b}{NOHD - a}\right]^{2} ; OD = - {}^{10}\log T$$
 target area ~ R²

Standard safety distance = 0.25 m for 10 sec exposure.

Laser treatment

In general: pulsed laser treatment will consist of separate trains of separate pulses



Continuous laser: 1 long train of 1 long pulse.

Example

LASER INPUT DATA	VALUE	
Wavelength	532	nm
(full) Divergence angle	0.0125	rad
Profile correlation factor	2	_
Waist(spot)diameter	400e-6	m
Distance laser to focus(waist)	0.3	m
Basic pulse: pulse power	1.000E+06	W
: pulse energy	5.000E-03	J
: duration	5.000E-09	S
: rep.freq.	1000	Hz
Train of pulses: power	5.000E+00	W
: energy	5.000E+01	J
: duration	1.000E+01	S
: rep.freq.	0.100	Hz
<pre>: pulse at end of train (yes/no=1/0)</pre>	0	-
: nr. of pulses per train	10000	-
Total treatment: power	5.000E+00	W
: energy	5.000E+01	J
: duration	1.000E+01	S
: nr. of trains in total	1	-
Safe distance from focus (for OD-glasses)	0.250	m

Example

```
OUTPUT of calculations:
For this wavelength the damage is:
EYE: photochemical and retina (thermal) and thermal
SKIN: thermal
            (H in J/m^2 - E=MPE in W/m^2)
RESULTS:
                formula
                                                                          NOHD [m] T-glasses OD-glass
                                                MPE [W/m^2] alpha [rad]
CALCULATION
 eye: pulse
                                                1.000E+06
                                                            3.133E-06
                                                                           127.9
                                                                                    3.837E-06
                                                                                                5.416
                H=5E-3.CE
 ..id. with CF1 .. N=
                                                1.000E+05
                                                            9.906E-07
                                                                           404.0
                                                                                    3.836E-07
                                                                                                6.416
                        10000, CF= 1.000E-01
                                                                           404.0
 ..id. with CF2 .. N=
                        10000, CF= 1.000E-01
                                                1.000E+05
                                                            9.906E-07
                                                                                    3.836E-07
                                                                                                6.416
                H=18.CE.t^0.75
                                                1.012E+01
                                                            4.461E-06
                                                                           89.99
                                                                                    7.769E-06
                                                                                                5.110
 eye: train
 eye: total
                                                                           89.99
                H=18.CE.t^0.75
                                                1.012E+01
                                                            4.461E-06
                                                                                    7.769E-06
                                                                                                5.110
                                                                           0.906
 skin: pulse
                H=200.CA
                                                4.000E+10
                                                            6.605E-04
                                                                                    1.700E-01
                                                                                                0.769
 ..id. with CF1 .. N=
                                                            2.012E-04
                                                                           2.286
                                                                                    1.584E-02
                                                                                                1.800
                        10000, CF= 1.000E-01
                                                4.000E+09
                        10000, CF= 1.000E-01
                                                                           2.286
                                                                                                1.800
 ..id. with CF2 .. N=
                                                4.000E+09
                                                            2.012E-04
                                                                                    1.584E-02
 skin: train
                H=1.1E4.CA.t^0.25
                                                1.956E+03
                                                            6.230E-05
                                                                           6.722
                                                                                    1.515E-03
                                                                                                2.819
 skin: total
                                                            6.230E-05
                                                                           6.722
                                                                                                2.819
                H=1.1E4.CA.t^{0.25}
                                                1.956E+03
                                                                                    1.515E-03
CF = correction factor, needed for thermal damage when >1 pulse in Tmin.
(Tmin is 1.8E-5 s for 532 nm).
```

the end